

LockSeq: Ensuring genomic safety in gene therapies through high-fidelity NGS validation


Felix Neumann

CEO & Co-Founder, Countagen AB

ICGT CONGRESS 2026 · Paris

May 27, 2026

Off-Target Confirmation Is Critical For Safety Assessment




- 01 CRISPR ideally generates **one intended clean edit** – at the same time, there are **hundreds of potential off-target sites** across the genome.
- 02 Current methods trade sensitivity for scale: **the more sites are surveyed the more dropouts***, leaving regulatory blind spots.
- 03 Off-targets detected late can **cost 6–50× more** to resolve than at pre-clinical stage.




* Dropout = The stochastic or systematic failure to capture a genomic locus during off-target validation.

LockSeq: High-Fidelity Molecular Lock Capture Chemistry

I. Target Capture



II. Library Preparation



One probe per target. No primer competition. No amplification bias.

Robust Variant Detection Across Platforms

0.1%

VAF sensitivity limit

= 1 in 1,000 alleles

$R^2 > 0.99$

ONT x Illumina x expected

Excellent correlation across platforms

>98%

Indel profiling accuracy

High confidence variant calling


ALSO VALIDATED

Cas9

Base editing

Prime editing

Transposon systems



OBSERVED VS. EXPECTED VAF · ILLUMINA × ONT

Acknowledgements

Minimal Target Dropout & Uniform Capture

0

LockSeq Dropouts

Across 101 designed sites

6

rhAmpSeq Dropouts

On same panel, same DNA

70%


Recovery Rate

Reduced Blind Spots


90%

Coverage Uniformity

Bias reduction



Site-retention through the panel pipeline



Recovery Rate Across 4 Panels


High Confidence Through Consensus Sequencing

Each captured site is tagged with a UMI. Consensus sequencing across reads grouped with UMI corrects stochastic and systematic errors – leaving only biologically true variants.

CAPTURED TARGETS



UMI GROUPED READS



CONSENSUS SEQUENCE



	WT	Edit
GLIS1	3.2	4.0
HCN1	4.6	16.8
HCP5B	4.4	5.3
IL18R1	6.6	6.8
ITS2	0.0	12.5

noisy, many false positives

	WT	Edit
GLIS1	0.0	0.0
HCN1	0.0	12.4
HCP5B	0.0	0.0
IL18R1	0.0	0.0
ITS2	0.0	0.0

clean, high-confidence variant calls

Two Ways To Start – Today

Full LockSeq Service

DNA-in → Report in 15 days. Custom panels up to 1,000+ sites.



- Bring DNA, walk away with data
- Sequencing platform of your choice

Strategic Partnerships

For teams thinking about integrating LockSeq in their own pipelines – long-term.

- Early-stage conversations open for facilities and CROs evaluating in-house deployment
- Karolinska CFG core facility already offers LockSeq for their customers – backend service offering
- Let's talk – we want to understand how can we support your off-target confirmation!

Trusted by



& more

Every edit.
Every site.
No blind spots.

Let's Connect!

Felix Neumann, CEO | felix@countagen.com

www.countagen.com

Stockholm, Sweden